#### CYBERGENE

# ChromoQuant® QF PCR kit Optima XY

## Detection of aneuploidy in chromosomes X and Y

## IVD kit for fast and accurate diagnosis of

- Aneuploidy in the sex chromosomes
- · Klinefelter syndrome (XXY)
- Turner syndrome (X0)
- Sex determination
- Deletion on the short arm of chromosome X

### Key advantages

- Single tube test with markers only for chromosomes X and Y
- 13 sex chromosome markers in total
- X Chr. markers on short and long arm of the chromosome
- 2 Y chromosome markers
- Turner X0 verifying TAF9B marker
- · Taq polymerase included, ready to use
- DNA source: Blood/Saliva/mouth swab/Amniotic fluid/CVS
- Detection of maternal contamination eliminates the risk of misdiagnosis
- Results are achieved within 6 hours enabling a "time to reply" of less than 24 hours.
- Used in prenatal testing or other types of human DNA tests
- The ChromoQuant® kits are validated for all sequencers from Life Technologies/
   ThermoFisher Scientific

#### **High flexibility**

The Optima XY kit can be used as a complement to the ChromoQuant AZF kit in order to diagnose Klinefelter samples when investigating male infertility. It can also be used for verification purpose in combination with the ChromoQuant Optima prenatal tests.

#### **High specificity**

The Optima XY kit holds 13 unique genetic markers.

#### **Interpretation of results**

GeneMarker and GeneMapper templates for easy interpretation and reporting are available as downloads.



#### CYBERGENE

# **Product specification**

CyberGene AB Banvaktsvägen 22 SE-171 48 Solna Sweden

Telephone: +46 8 608 23 90

Sales support: sales@cybergene.se

Technical support: chromoquant@cybergene.se

www.cybergene.se

Intended Use	In vitro diagnostics for diagnosis of chromosome X and Y aneuploidy			
No. of markers in Optima XY	-	<b>Y</b> 2	<b>XY</b> 4	Turner 1
Complies with Best Practice guidelines for QF-PCR	Yes			
CE-labelled for IVD use	Yes			
Taq polymerase	Included			
Detection format	Capillary Electrophoresis with Genetic Analyser			
Validated Genetic Analysers	ABI 310, 3100, 3130, 3730, 3500, SeqStudio			
Part no. and kit size	514	4.60	4-26 (	(26 tests)



ChromoQuant® has been thoroughly validated. ChromoQuant® was clinically introduced in early 2004 and is used world wide.

ChromoQuant $^{\otimes}$  is CE marked in accordance with the IVD Directive 98/79/EC and produced by CyberGene AB under quality system ISO 13485.

#### **About CyberGene AB**

CyberGene AB is a Swedish company, active in the MedTech field by developing, manufacturing and selling In Vitro diagnostic products. ChromoQuant is a registered trademark of CyberGene AB.

